


APPARATUS 1 & 2


QLA manufactures a full line of standard 1000mL nominal capacity **Dissolution Vessels**.

These vessels are made of high quality borosilicate glass. 2000mL and 500mL nominal size vessels are also available. All glass vessels have graduations and serial numbers.

All QLA dissolution vessels are USP compliant and meet machine manufacturer specifications



**GLA900-01
GLA900-HR
GLA900-DK**




GLA900-AD

GLA900-AV

GLA900-AH

GLA500-HR

GLA500-01



GLA900-VK

**PLA900-01
PLA900-HR
PLA900-DK**

**PLA900-AV
PLA900-AH
PLA900-AD**

QLA also carries a clear plastic vessel which fits all major testing machines as well as the Hanson Mediamate machines. Plastic vessels can also be tinted amber for light sensitive applications.

APPARATUS 1 & 2




QLA P/N	Product Description	VanKel P/N	Hanson P/N	Distek P/N
GLA900-01	1000mL Clear Glass Vessel	12-5000		
GLA900-AV	1000mL Amber Glass Vessel	12-5010		
GLA900-VK	1000mL Nominal "V Series" Glass Dissolution Vessel	12-5035		
GLA02K-01	2000mL Clear Glass Vessel	12-5070		
GLA500-01	500mL Clear Glass Vessel	12-5060		
PLA900-01	1000mL Clear Plastic Vessel	12-5200		
PLA900-AV	1000mL Amber Plastic Vessel	New QLA innovated product used in VanKel applications		
VK600A-CR	Plastic Centering Ring Assembly for VK600/6000 Series	12-6000		
GLA900-HR	1000mL Clear Glass Vessel		72-600-510	
GLA900-AH	1000mL Amber Glass Vessel		72-600-556	
GLA500-HR	500mL Clear Glass Vessel	New QLA innovated product used in Hanson applications		
PLA900-AH	1000mL Amber Plastic Vessel	New QLA innovated product used in Hanson applications		
PLA900-HR	1000mL Clear Plastic Vessel		91-030-003	
GLA200-DK	200mL Clear Glass Vessel - No Ring	New QLA innovated product used in Distek applications		
GLA900-DK	1000mL Clear Glass Vessel	QLA offers this item without centering ring for cost savings		
GLA9CR-DK	1000mL Clear Glass Vessel	with Centering Ring		
GLA900-AD	1000mL Amber Glass Vessel			3010-0093
GLA9CR-AD	900mL Amber Glass Vessel with Centering Ring	New QLA innovated product used in Distek applications		
GLA02K-DK	2000mL Clear Glass Vessel			3010-0285
GLA500-DK	500mL Clear Glass Vessel	New QLA innovated product used in Distek applications		
PLA900-DK	1000mL Clear Plastic Vessel			3350-0095
PLA900-AD	1000mL Amber Plastic Vessel	New QLA innovated product used in Distek applications		
RING01-DK	Plastic Centering Ring (Acculign)			3350-0096

PEAK VESSELS

GLA9PK-01	1000mL PEAK Glass Vessel	12-5500
-----------	--------------------------	---------

QLA supplies patented **PEAK Vessels** which are designed to prevent cone formation at the vessel base.

Clear PEAK Vessels are always in stock. These vessels can be Amber or PTFE coated as special orders



GLA9PK-01