







800 Series MICROLITER Syringes

5 µL - 250 µL

- Removable needles (RN) or cemented needles (N)
- Reinforced plungers
- Plunger assemblies are replaceable
- Plunger stop prevents plunger blowout
- For use with liquids
- Syringes are autoclavable when disassembled except for cemented needle (N) syringes

Ordering Information

		Volume	5 µL	10 µL	25 µL	50 µL	100 µL	250 µL
		Model	85	801	802	805	810	825
		Gauge **	26s	26s	22s	22s	22s	22s
N (Cemented Needle)			84850	84852	84854	84856	84858	
RN (Removable Needle)			84851	84853	84855	84857	84859	84861
RN (Removable Needle) Syringe (Needle not included)			7641-01	7642-01	7643-01	7644-01	7645-01	7646-01
Reproducibility (Chaney) Adapter (Syringe not included)			32146	32146	32146	32146	32146	32146



Syringe Calibration Service


Add the prefix "CAL" to the beginning of the syringe's part number.

Typical Operating Parameters

- Maximum temperature for Cemented Needle 50 °C (122 °F) and Removable Needle 115°C (239°F).
- Maximum pressure rating is 2000 psig except for 250 µL which is 1000 psig.

Replacement Parts

Description	Volume	5 µL	10 µL	25 µL	50 µL	100 µL	250 µL
	Model	85	801	802	805	810	825
1. Barrel/Plunger Assembly N		32150	32151	32152	32153	32154	32169
Barrel/Plunger Assembly RN		32134	32129	32117	32120	32123	32126
2. Holder Assembly		32135	32135	32135	32135	32135	32135
3. RN Nut		30902	30902	30902	30902	30902	30902
4. RN Needle, see page 31							









900 Series MICROLITER Syringes

5 µL - 10 µL

- Removable needles (RN) or cemented needles (N)
- Reinforced plungers
- For use with liquids
- Syringes are autoclavable when disassembled except for cemented needle (N) syringes

Ordering Information

		Volume	5 µL	10 µL
		Model	95	901
		Gauge **	26s	26s
N (Cemented Needle)			87920	80360
RN (Removable Needle)			87925	80370
RN (Removable Needle) Syringe (Needle not included)			7647-01	7648-01
Reproducibility (Chaney) Adapter (Syringe not included)			32148	32148



Syringe Calibration Service

Add the prefix "CAL" to the beginning of the syringe's part number.

Typical Operating Parameters

- Maximum temperature for Cemented Needle 50 °C (122 °F) and Removable Needle 115 °C (239 °F).
- Maximum pressure rating is 2000 psig.

Replacement Parts

Description	Volume	5 µL	10 µL
	Model	95	901
1. RN Nut		30902	30902
2. RN Needle, see page 31			

